
FREND Free T4

REF FRT4AP 020

- English
- German
- Italian
- Spanish
- French
- Portuguese

FREND Free T4

Quantitative Assay for Free Thyroxine

REF FRT4AP 020

IVD For *in vitro* diagnostic use only

1. Intended use

The FREND Free T4 Test System is a rapid indirect competitive fluorescence immunoassay for the quantitative determination of free thyroxine (FT4) in human serum and lithium heparinized plasma specimens using the FREND™ System. Measurements of FT4 are used in the diagnosis of thyroid disorders. The FREND Free T4 Test is intended for use in clinical laboratories. For *in vitro* diagnostic use only. For professional use only.

2. Summary and explanation of test

Human Thyroid Stimulating Hormone (hTSH) or thyrotropin stimulates the secretion of thyroxine (T4) and triiodothyronine (T3) by the thyroid gland.¹

From the moment it is secreted into the bloodstream, thyroxine or tetra iodothyronine (T4), produced by the thyroid gland, is predominantly (>99%) bound to the carrier proteins TBG (Thyroxine Binding Globulin), TBPA (Thyroxine Binding Pre Albumin) and albumin. The fraction that remains free (FT4) is considered the active part of the hormone.² The mechanisms regulating thyroid function have a direct effect on the concentration of this free fraction, which explains why it is relatively independent of the concentration of carrier proteins.^{3,4} In patients with hyperthyroidism, the FT4 concentration increases, whereas in patients with hypothyroidism it generally decreases. Patients on levothyroxine hormone replacement therapy (LT4)

may have an elevation of FT4, although clinically they are euthyroid.

T3 and T4 are known to have diverse functions in regulating basal metabolic rate, bone growth, neuronal development, and sexual maturation.⁵ Underproduction of T3 and/or T4 can result in hypothyroidism, while overproduction of these hormones can result in hyperthyroidism.⁶ Because of a negative feedback mechanism, elevation of T3 and T4 suppress the production of TSH.⁷ TSH itself is stimulated by thyrotropin-releasing hormone (TRH), a tripeptide produced in the hypothalamus.⁸

Primary hypothyroidism occurs when TSH levels are elevated while T3 and/or T4 are underproduced.^{9,10} Secondary or tertiary hypothyroidism can occur due to an abnormal response of TSH to TRH, while central hypothyroidism occurs from pituitary dysfunction.^{11,12} Primary hyperthyroidism is marked by low levels of TSH and high levels of T3 and/or T4.⁹ Anomalies to these types of classification exist, but TSH testing can (with the aid of other thyroid tests) help a clinician determine the presence of thyroid dysfunction.¹³⁻¹⁵

3. Principle of the assay

A 70µl specimen is added to pretreatment tube containing gold-t4 antibodies and incubated in the FREND™ AP System automatically, while free thyroxines react with gold particles bound to antibodies against free thyroxine.

After five minutes, the AP device transfers a 35 µL of sample to FREND Free T4 cartridge and incubated in AP system automatically. The cartridge is then placed into the FREND™ System, which is programmed to begin analysis once the sample has reacted with the reagents. The reaction and analysis time is ample, so a lower ratio of fluorescence correlates with a higher FT4 concentration.

The FREND™ System is a bench-top fluorescence reader containing a touchscreen user interface. The System has a slot that accepts the FREND Free T4 test cartridge (which contains the reagents and sample), and is programmed to analyze the test when the sample has fully reacted with the on-board cartridge reagents. Results of the test are displayed on the screen and can be printed on an optional printer.

4. Material provided

Q'ty	Contents	Catalogue number
20	FREND Free T4 Cartridges	FRT4AP 020
20	FREND Free T4 Pretreatment tubes	
30	Disposable pipette tip(s)	
01	FREND Free T4 Code chip	
01	FREND Free T4 Package insert	

One cartridge contains:

T3-BSA
Anti-T3 antibody
Fluorescent particles

One pretreatment tube contains:

Anti-T4 antibody
Gold nano-particle

5. Materials required but not provided

- The FREND™ System
- The FREND™ AP
- Micro-pipette, or any pipette, capable of delivering 35 and 70 µL
- Personal protective equipment and biohazard waste disposal containers

6. Warning and Precautions

Caution: Federal law restricts this device to sale by or on the order of a physician.

- The FREND Free T4 cartridges are intended for in vitro diagnostic use only.
- The FREND Free T4 cartridges are only to be used on the FREND™ System.
- The FREND Free T4 cartridges are disposable, single use devices. Do not reuse them under any circumstances.
- Allow sealed cartridges to come to room temperature for 15-30 minutes prior to use.
- Cartridges and pretreatment tubes should not be frozen.

- Assure the humidity in the laboratory is in the 10-80% range and that the room temperature is in the range of 64–77°F (18–25°C) when tests run.
- Avoid cross-contamination between samples by using a new pipette tip for each new specimen.
- Avoid high humidity, direct sunlight or heat in the area used for cartridge storage.
- Inaccurate results are possible if the sample used is contaminated in any way.
- Using specimens containing clotted fibrin could result in erroneous results
- Over or under loading the cartridge with sample may result in inaccurate results.
- Do not use the cartridges beyond the expiration date on the pouch.
- Do not use the cartridge and pretreatment if the pouch is damaged or the seal is broken.
- Perform testing as specified in the Package Insert Sheet and User Manual.
- Keep the cartridge and pretreatment tube sealed in the pouch until ready for use.
- Use the cartridge immediately after opening the pouch.
- Handle specimens in accordance with the OSHA Standard on Bloodborne Pathogens.
- Human specimens are not used in the preparation of this product, however, since human specimens will be used for samples and other quality control products in the laboratory may be derived from human materials, Use Universal Precautions when handling all specimens and controls. Wear disposable gloves when handling the cartridges and the samples.
- Wash hands thoroughly and often after handling reagent cartridges or samples.
- Do not ingest the silica gel packet found in the cartridge pouch.

7. Storage and Stability

All unopened materials are stable until the expiration date on the label when stored at the specified temperature. Reagent stability has been demonstrated for 18 months from the date of manufacture.

The expiration date is clearly indicated on the product box and the cartridges

Materials	Catalogue number
Refrigerator temperature storage (2–8 °C)	
FREND Free T4 Cartridges	FRT4AP 020
FREND Free T4 Pretreatment tubes	None
Room temperature storage (18–25 °C)	
Pipette tips	None

8. Specimen collection and handling

Serum or lithium heparinized plasma is required for the assay. No special patient preparation is necessary. Collect the appropriate venous blood sample aseptically. For serum, allow the sample to clot for 30 minutes at room temperature. For lithium heparin, centrifuge after collection. Centrifuge the sample for 10 minutes at 3,000 rpm within 2 hours of collection and immediately separate the serum or plasma from the packed cells.

Separated Samples may be stored at 2-8 °C for up to 1 week prior to analysis. If the analysis is scheduled to be done more than 1 week after collection, the sample should be stored frozen at -20 °C or below for future use.

Repeated freeze-thaw cycles should be avoided. Turbid serum samples or samples containing particulate matter such as fibrin clots or visible strands should be re-centrifuged before being tested. Prior to assay, slowly bring frozen samples to room temperature and mix gently but thoroughly before testing.

The sample required for the incubation step is 70 µL. The sample required for running the test on the FREND Free T4 cartridge is 35 µL.

For optimal results, avoid grossly hemolytic, lipemic or turbid specimens. Specimens should be free of aggregated fibrin, red blood cells, or other particulate matter. When pipetting into the FREND Free T4 cartridge sample inlet, ensure that bubbles in the sample are avoided. Bubbles may restrict flow and result in an incomplete or erroneous test result.

9. Procedure

Calibration

There is no need for calibration to be performed by the end user. All calibration statistics and information have been electronically stored on the FREND Free T4 Code chip included in each box of FREND Free T4 cartridges. The FREND Free T4 Code chip is specific for each manufactured lot of FREND Free T4 cartridges.

Always run external quality control samples to verify that the FT4 results obtained on the FREND™ System meet the laboratory criteria for acceptability for each lot of FREND Free T4 cartridges.

Code chip installation

Please refer to the FREND™ System User manual for more detailed instructions relative to the Code chip installation. Abbreviated instructions follow here:

- 1) Insert the FREND™ System electrical cord into an appropriate outlet.
- 2) Insert the Code chip into the Code chip slot at the rear of the FREND™ System following the arrows.
- 3) Press the **'Setup'** button on the **'Main'** screen.
- 4) Press the **'Code chip'** button on the **'Setup'** screen.
- 5) The information embedded on the FREND Free T4 Code chip is automatically saved on the FREND™ System.
- 6) When the Code chip installation is completed, press the **'OK'** button to go to the **'Setup'** screen.
- 7) Press the **'Item'** button on the **'Setup'** screen.
- 8) Check the FREND Free T4 cartridge lot number and the installation date of the Code chip.
- 9) Press the **'Home'** button to go to the **'Main'** screen to begin running external quality control and patient samples.

Quality control

• FREND™ System QC cartridges

FREND QC Cartridge contains multiple controls to check optic part of the system. By testing QC Cartridge, part of analytical components of the system of (1) laser power, (2) alignment, and (3) mechanical integrity are confirmed.

For each day of patient testing, perform QC Cartridge testing. Refer to the quality control procedures section in the User Manual of FREND™ System. In brief, perform QC cartridge testing for the following conditions

- 1) Upon initial setup of the system
- 2) Each day of patient testing,
- 3) When the system has been transported or moved,
- 4) Whenever there is uncertainty about the performance of the system,
- 5) Whenever required by your laboratory's quality control requirements.

• **Internal procedural controls**

The FREND Free T4 test cartridge contains built-in control features. Fluorescence signal in the Reference Zone of each cartridge shows: (1) that enough sample volume is added, (2) that proper flow is obtained, and (3) that the antibody is reactive. If this Reference Zone signal is missing or lower than the threshold, the FREND™ System considers it as an incorrect or failed test, and produces an error message instead of a test result. In addition, with each cartridge run, the system monitors, in part, for (1) flow of sample, (2) speed of sample flow, (3) shelf-life of cartridge components, (4) function of internal barcode scanner, and (5) function of scanner's mechanical components.

• **External quality control testing**

Commercially available controls that contain FT4 as a measured analyte are available from a variety of manufacturers. It is recommended that a minimum of two (2) levels of controls be run at least once per month or once for each new lot, whichever comes earlier. However, Controls should be run according to the local requirements for each laboratory. Each laboratory should establish its own criteria based on the following parameters.

- 1) Each new lot,
- 2) Each new shipment (even if from the same lot previously received),
- 3) Each new operator (an individual who has not run the tests for at least two weeks),
- 4) Monthly, as a continued check on storage conditions,
- 5) Whenever problems (storage, operator, or other) are identified,
- 6) Or other times as required by your laboratory's standard QC procedures.

Individual laboratory policy will dictate exactly which control materials and lot numbers should be run, the frequency with which controls are to be tested criteria for acceptance if the results and required corrective action to be taken if results do not meet laboratory criteria. If any external quality control sample values are out of the acceptable range, it will be necessary to investigate the problem before reporting patient results to assure there is not an instrument of software malfunction. Do not assay patient samples on the FREND™ System using FREND Free T4 if quality control results do not give expected

values. Refer to your laboratory policies on how to determine acceptability of external control materials results. Every laboratory must follow the standardized procedures acceptable to the regulatory agencies to whom the laboratory is responsible.

Specimen processing

• Preparation

Remove sufficient FREND Free T4 cartridges and pretreatment tubes from the refrigerator to test the number of patient samples and required external quality control materials. Allow the tubes and the sealed pouches containing the cartridges to come to room temperature for 15-30 minutes prior to the start of the testing sequence.

If using refrigerated patient samples, remove those from the refrigerator and allow them to come to room temperature prior to testing. If frozen samples will be utilized, be sure these are removed from the freezer, thawed naturally and then mixed gently but thoroughly prior to testing. Testing should not begin on previously frozen samples until they have reached room temperature.

There are no other reagents or sample preparations necessary.

• Assay procedure

Note: *When processing samples, rinse the pipette several times with the sample and dispense 70 µL into the pretreatment tube.*

- 1) Prepare the FREND Free T4 cartridge, pretreatment tube and specimen at room temperature (18–25°C). Open the pouch and place the FREND Free T4 cartridge into the cartridge tray of the AP device. Press “**NEXT**” to close the cartridge tray and open the pretreatment tube tray.
- 2) Transfer 70 µL of specimen to the pretreatment tube.

Caution: *Once the sample is added to the pretreatment tube, do not invert the tube.*

- 3) Insert the pretreatment tube into the tube hole in the FREND™ AP pretreatment tube tray. Refer to the FREND™ AP User manual for complete operating instructions.
- 4) Press the “**NEXT**” button. The pretreatment tray will close and the first incubation step (5 minutes) will begin.

- 5) After the first incubation is completed, 35µL of mixed sample will be loaded onto the cartridge and the second incubation step (2 minute) will begin.
- 6) When both incubation steps are completed, the cartridge tray will open and the cartridge will be ready to be inserted into the FREND™ System.
- 7) Press the '**Test**' button on the '**Main**' screen of the FREND™ System.
- 8) The system moves to the Patient ID screen automatically.
- 9) Type the Patient ID and press the '**Enter**' button to begin the test.
- 10) Insert the cartridge into the cartridge slot using the cartridge arrow as a guide.

Caution: Check the direction of the cartridge before insertion and assure the insertion is complete.

- 11) When the reaction in the cartridge is completed, the FREND™ System will automatically begin the reading process.
- 12) When the measurements are completed, the cartridge will automatically be expelled and the results displayed.

Caution: Do not disconnect power cord or shut off power on the FREND™ System while a cartridge is in the reading chamber. This may cause a system error.

- 13) If the FREND™ System is connected to the optional printer, press the 'Print' button and the results will be output on the printer paper.
- 14) For more detailed instructions, please refer to the FREND™ System User Manual.

10. Specimen Dilution Procedures

Samples cannot be diluted for FT4 determinations. Samples which read >6.00 ng/dL should be reported as such.

11. Calculation of results

The FREND™ System performs all sample and reagent handling operations automatically within the cartridge once the sample has been manually loaded to the sample inlet in the cartridge and the cartridge placed into the FREND™ System. The rate of fluorescence produced by the reaction is read at various intervals during the analysis process, blank readings are subtracted after which the net rate is automatically converted to FT4 concentration in ng/dL based upon information stored on the Free T4 Code chip. This result is then displayed on the screen and can be sent to the optional printer. It is also stored in memory on the FREND™ System.

Screen displayed for various concentration scenarios

Displayed result	Description
	FT4 Concentration Less than 0.40 ng/dL
	FT4 Concentration Not less than 0.40 ng/dL And not higher than 6.00 ng/dL
	FT4 Concentration Higher than 6.00 ng/dL

12. Limitations of the procedure

- 1) When used for diagnostic purposes, the results obtained from this assay should be used in conjunction with other data (e.g, symptoms, results of other tests, clinical impressions, medical history, therapy, etc.)
- 2) The FRENDS™ System paired with a FRENDS Free T4 cartridge, is programmed to report 6.00 ng/dL as the highest concentration of FT4 measurable without dilution. The lowest measurable concentration is 0.40 ng/dL – the assay limit of quantitation.

- 3) Specimens from patients with heterophilic antibodies, such as anti-mouse (HAMA), anti-goat (HAGA), or anti-rabbit (HARA) antibodies, may show falsely elevated or depressed values or may result in the error message **“Incomplete Test”**. Patients routinely exposed to animals or animal serum products can be prone to these types of heterophilic interferences. If the FT4 level is inconsistent with clinical evidence, additional FT4 or other thyroid testing using a different method is suggested to confirm the results.
- 4) Certain medications may interfere with assay performance. All results should be interpreted with respect to the clinical picture of the patient.
- 5) Although hemolysis has an insignificant effect on the assay, hemolyzed samples may indicate mistreatment of a specimen prior to assay and results should be interpreted with caution.
- 6) Lipemia has an insignificant effect on the assay except in the case of gross lipemia where interference with the lateral flow of the sample in the cartridge may occur.
- 7) The concentration of FT4 in a given sample determined using assays from different manufacturers can vary due to differences in assay methods, calibration, and antibody specificity.
- 8) Please refer to the Specimen Collection and Handling, Warnings and Precautions, Storage and Stability, and Procedural Notes sections in this insert sheet.
- 9) FREND Free T4 has not been validated in point-of-care settings.
- 10) Performance of this assay has not been established with neonatal specimens or specimens from pregnant patients.
- 11) FREND Free T4 is to be used in licensed clinical laboratories with trained technologists.

13. Expected values

As with every clinical diagnostic test, a reference interval corresponding to the characteristics of the population being tested should be determined by each laboratory. Historically, it has been shown that there are no race- or gender-specific differences in the reference interval for FT4, so a single adult reference interval is reasonable and justified. Serum samples from a total of 196 normal, apparently healthy adult individuals were assayed on 3 lots of the FREND Free T4 assay using a single FREND™ System. The reference interval was determined according to CLSI guideline C28-A3, was found to be 0.83-1.60 ng/dL. As in all in vitro diagnostic testing, a Free T4 result generated using the FREND Free T4 on the FREND™ System should be interpreted in the light of other clinical findings and diagnostic procedures. Any FT4 results not correlating with the clinical condition should be repeated and other testing performed to clarify the situation.

14. Performance characteristics

Precision


A precision study for FREND Free T4 was performed as outlined in CLSI EP5-A3. Four samples were assayed in two replicates, twice per day over a period of 20 days.

Sample ID	Mean Conc. (ng/dL)	Within-run		Between-run		Between-day		Within-laboratory	
		SD	CV(%)	SD	CV(%)	SD	CV(%)	SD	CV(%)
1	0.807	0.066	8.2	0.013	1.6	0.015	1.9	0.069	8.6
2	1.812	0.140	7.7	0.033	1.8	0.009	0.5	0.145	8.0
3	2.513	0.116	4.6	0.096	3.8	0.016	0.6	0.151	6.0
4	4.060	0.178	4.4	0.083	2.1	0.065	1.6	0.207	5.1

Linearity

The dilution linearity study as outlined in CLSI EP6-A was performed by dilution a high concentration FT4 specimen with FT4 depleted serum. Linearity was demonstrated from 0 to 6 ng/dL.


FREND Free T4 - Linearity

The analytical measurement range of FREND Free T4 is from 0.40 ng/dL to 6.00 ng/dL.

Method comparison

FREND Free T4 was compared to the predicate device using guidelines outlined in CLSI document EP09-A3. Samples (n=120) were measured in duplicate on both systems. Linear regression analysis demonstrated a correlation coefficient (R^2) of ≥ 0.98 .

Sensitivity

The limit of blank (LoB) and limit of detection (LoD) was determined using guidelines found in CLSI document EP17-A. LoB was determined from 60 replicate measurements using a negative standard solution. LoD was determined using 12 replicates measurements of five low patient samples.

LoB	LoD	LoQ
0.29 ng/dL	0.36 ng/dL	0.36 ng/dL

Interference

The interference study was performed as recommended in the CLSI EP7-A2 protocol using three concentrations of FT4. Recovery within 90% to 110% of the expected FT4 was considered as lack of interference. No interference by the substances below was found.

No.	Potential interfering	Concentration
1	Hemoglobin	500 mg/dL
2	Bilirubin	20 mg/dL
3	Triglycerides	3 g/dL
4	Total protein	12 g/dL

Cross-reactivity

The following substances were evaluated for potential cross-reactivity with the FREND Free T4 at two concentrations. Testing was done according to the CLSI protocol EP7-A2. No significant cross-reactivity was found.

No.	Potential cross reactants	Concentration(ng/dL)
1	Levothyroxine, T4	1,000
2	Diiodothyronin, T2	5,000
3	Tetraiodothyroacetic Acid	10,000
4	Triiodothyroacetic Acid	1,000
5	Triiodothyropropionic Acid	5,000
6	Diiodothyrosine, DIT	1,000,000
7	L-Triiodothyronine, T3	1,000
8	Monoiodotyrosine	1,000,000
9	Reverse T3	10,000

15. References


- 1) Pierce JG. The subunits of pituitary thyrotropin—their relationship to other glycoprotein hormones. *Endocrinology* 1971, 89, 1331-1344.
- 2) Sapin R, D'herbomez M. Dosage des hormones thyroïdiennes: thyroxine (T4) et triiodothyronine (T3). In *La thyroïde*, LeClere J, Orgiazzi J, Rousset B, Schlienger Elsevier, ed. 2001, 268-274.
- 3) Pearce C.J. Byfield PGH. Free thyroid hormone assays and thyroid function. *Annual Clinical Biochemistry*. 1986, 23, 230-237.
- 4) Ekins R. Measurement of Free Hormones in Blood - *Endocrine Reviews*, Vol. 11, No. 1, 1990, 5-46.
- 5) Smith BR, Pyle GA, Peterson VB, Hall R. Interaction of thyrotropin with the human thyrotropin receptor. *J. Endocrinol.* 1977, 75, 391-400.
- 6) Sterling K, Lazarus JH. The thyroid and its control. *Annu. Rev. Physiol.* 1977, 39, 349-371.
- 7) Patel YC, Alford FP, Burger HG. The 24-hour plasma thyrotrophin profile. *Clin. Sci.* 1972, 43, 71-77.
- 8) Morley JE. Neuroendocrine control of thyrotropin secretion. *Endocr. Rev.* 1981, 2, 396-436.
- 9) Wehmann RE, Rubenstein HA, Pugeat MM, Nisula BC. Extended Clinical Utility of a Sensitive and Reliable Radioimmunoassay of Thyroid-Stimulating Hormone. *South. Med. J.* 1983, 76, 969-976.
- 10) Burger HG, Patel YC. The Value of Serum Thyrotropin Measurement in the Diagnosis and Management of Hypothyroidism. *Med. J. Aust.* 1972, 2, 293-297.
- 11) Petersen VB, McGregor AM, Belchetz PE, Elkeles RS, Hall R. The Secretion of Thyrotropin with Impaired Biological Activity in Patients with Hypothalamic-Pituitary Disease. *Clin. Endocrinol.* 1978, 8, 397-402.
- 12) Beck-Peccoz P, Amr S, Menezes-Ferreira MM, Faglia G, Weintraub B. Decreased Receptor Binding of Biologically Inactive Thyrotropin in Central Hypothyroidism. *N. Engl. J. Med.* 1985, 312, 1085-1090.
- 13) Hay ID, Bayer MF, Kaplan MM, Klee GG, Larson PR, Spencer CA. American Thyroid Association Assessment of Current Free Thyroid Hormone and Thyrotropin Measurements and Guidelines for Future Clinical Assays. *Clin. Chem.* 1991, 37, 2002-2008.
- 14) Demers LM, Spencer CA. Laboratory Medicine Practice Guidelines: Laboratory Support for the Diagnosis and Monitoring of Thyroid Disease. *The National Academy of Clinical Biochemistry*, 1996.
- 15) Beastall GH, Beckett GJ, Franklyn J; Frasier WD, Hickey J, John R, Kendall-Taylor P, Nevens B, Vanderpump M. UK Guidelines for the Use of Thyroid Function Tests. *The Association for Clinical Biochemistry*, 2006.

Glossary of Symbols

	Caution, warning, Consult accompanying documents
	Catalogue number/Reference number
 www.nanoentek.com/eifu.php	Consult Instructions for Use An electronic instructions for use (eIFU) indicator (website address) may accompany the symbol when used to indicate an instruction to consult an eIFU.
	Lot number/Batch number
	Use by YYYY-MM-DD or YYYY-MM
	Manufacturer
	CE marking
	<i>In vitro</i> diagnostic medical device
	Temperature limitation
	Contains sufficient for <n> tests
	Do not reuse
	Do not use if package is damaged
	For prescription use only CAUTION: Federal (U.S.) law restricts this device to sale by or on order of a physician.
	US Corporation
	Patient ID
	Result
 Sample Drop	Sample Drop
	Authorized representative in the European Community

UK Representative	Authorized representative in United Kingdom
CH REP	Authorized representative in Switzerland
BRH	Authorized representative in Brazil

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15. Referencias

- 1) Pierce JG. The subunits of pituitary thyrotropin—their relationship to other glycoprotein hormones. *Endocrinology* 1971, 89, 1331-1344.
- 2) Sapin R, D'herbomez M. Dosage des hormones thyroïdiennes: thyroxine (T4) et triiodothyronine (T3). In La thyroïde, LeClere J, Orgiazzi J, Rousset B, Schlienger Elsevier, ed. 2001, 268-274.
- 3) Pearce C.J. Byfield PGH. Free thyroid hormone assays and thyroid function. *Annual Clinical Biochemistry*. 1986, 23, 230-237.
- 4) Ekins R. Measurement of Free Hormones in Blood - *Endocrine Reviews*, Vol. 11, No. 1, 1990, 5-46.
- 5) Smith BR, Pyle GA, Peterson VB, Hall R. Interaction of thyrotropin with the human thyrotropin receptor. *J. Endocrinol.* 1977, 75, 391-400.
- 6) Sterling K, Lazarus JH. The thyroid and its control. *Annu. Rev. Physiol.* 1977, 39, 349-371.
- 7) Patel YC, Alford FP, Burger HG. The 24-hour plasma thyrotrophin profile. *Clin. Sci.* 1972, 43, 71-77.
- 8) Morley JE. Neuroendocrine control of thyrotropin secretion. *Endocr. Rev.* 1981, 2, 396-436.
- 9) Wehmann RE, Rubenstein HA, Pugeat MM, Nisula BC. Extended Clinical Utility of a Sensitive and Reliable Radioimmunoassay of Thyroid-Stimulating Hormone. *South. Med. J.* 1983, 76, 969-976.
- 10) Burger HG, Patel YC. The Value of Serum Thyrotropin Measurement in the Diagnosis and Management of Hypothyroidism. *Med. J. Aust.* 1972, 2, 293-297.
- 11) Petersen VB, McGregor AM, Belchetz PE, Elkeles RS, Hall R. The Secretion of Thyrotropin with Impaired Biological Activity in Patients with Hypothalamic-Pituitary Disease. *Clin. Endocrinol.* 1978, 8, 397-402.
- 12) Beck-Peccoz P, Amr S, Menezes-Ferreira MM, Faglia G, Weintraub B. Decreased Receptor Binding of Biologically Inactive Thyrotropin in Central Hypothyroidism. *N. Engl. J. Med.* 1985, 312, 1085-1090.
- 13) Hay ID, Bayer MF, Kaplan MM, Klee GG, Larson PR, Spencer CA. American Thyroid Association Assessment of Current Free Thyroid Hormone and Thyrotropin Measurements and Guidelines for Future Clinical Assays. *Clin. Chem.* 1991, 37, 2002-2008.
- 14) Demers LM, Spencer CA. Laboratory Medicine Practice Guidelines: Laboratory Support for the Diagnosis and Monitoring of Thyroid Disease. *The National Academy of Clinical Biochemistry*, 1996.
- 15) Beastall GH, Beckett GJ, Franklyn J; Frasier WD, Hickey J, John R, Kendall-Taylor P, Nevens B, Vanderpump M. UK Guidelines for the Use of Thyroid Function Tests. *The Association for Clinical Biochemistry*, 2006.

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15. Références

- 1) Pierce JG. The subunits of pituitary thyrotropin—their relationship to other glycoprotein hormones. *Endocrinology* 1971, 89, 1331-1344.
- 2) Sapin R, D'herbomez M. Dosage des hormones thyroïdiennes: thyroxine (T4) et triiodothyronine (T3). In *La thyroïde*, LeClere J, Orgiazzi J, Rousset B, Schlienger Elsevier, ed. 2001, 268-274.
- 3) Pearce C.J. Byfield PGH. Free thyroid hormone assays and thyroid function. *Annual Clinical Biochemistry*. 1986, 23, 230-237.
- 4) Ekins R. Measurement of Free Hormones in Blood - *Endocrine Reviews*, Vol. 11, No. 1, 1990, 5-46.
- 5) Smith BR, Pyle GA, Peterson VB, Hall R. Interaction of thyrotropin with the human thyrotropin receptor. *J. Endocrinol.* 1977, 75, 391-400.
- 6) Sterling K, Lazarus JH. The thyroid and its control. *Annu. Rev. Physiol.* 1977, 39, 349-371.
- 7) Patel YC, Alford FP, Burger HG. The 24-hour plasma thyrotrophin profile. *Clin. Sci.* 1972, 43, 71-77.
- 8) Morley JE. Neuroendocrine control of thyrotropin secretion. *Endocr. Rev.* 1981, 2, 396-436.
- 9) Wehmann RE, Rubenstein HA, Pugeat MM, Nisula BC. Extended Clinical Utility of a Sensitive and Reliable Radioimmunoassay of Thyroid-Stimulating Hormone. *South. Med. J.* 1983, 76, 969-976.
- 10) Burger HG, Patel YC. The Value of Serum Thyrotropin Measurement in the Diagnosis and Management of Hypothyroidism. *Med. J. Aust.* 1972, 2, 293-297.
- 11) Petersen VB, McGregor AM, Belchetz PE, Elkeles RS, Hall R. The Secretion of Thyrotropin with Impaired Biological Activity in Patients with Hypothalamic-Pituitary Disease. *Clin. Endocrinol.* 1978, 8, 397-402.
- 12) Beck-Peccoz P, Amr S, Menezes-Ferreira MM, Faglia G, Weintraub B. Decreased Receptor Binding of Biologically Inactive Thyrotropin in Central Hypothyroidism. *N. Engl. J. Med.* 1985, 312, 1085-1090.
- 13) Hay ID, Bayer MF, Kaplan MM, Klee GG, Larson PR, Spencer CA. American Thyroid Association Assessment of Current Free Thyroid Hormone and Thyrotropin Measurements and Guidelines for Future Clinical Assays. *Clin. Chem.* 1991, 37, 2002-2008.
- 14) Demers LM, Spencer CA. Laboratory Medicine Practice Guidelines: Laboratory Support for the Diagnosis and Monitoring of Thyroid Disease. The National Academy of Clinical Biochemistry, 1996.
- 15) Beastall GH, Beckett GJ, Franklyn J, Frasier WD, Hickey J, John R, Kendall-Taylor P, Nevens B, Vanderpump M. UK Guidelines for the Use of Thyroid Function Tests. The Association for Clinical Biochemistry, 2006.

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NanoEntek, Inc.
MT Promedt Consulting GmbH