
FREND Vitamin D

 FRVDAP 020

- English
- German
- Italian
- Spanish
- French
- Portuguese

NESPI-FVAP-001ML(V.0.0)

FREND Vitamin D

Quantitative Assay for Total 25-Hydroxy Vitamin D

FRVDAP 020

*For in vitro diagnostic use only*

1. Intended use

The FREND Vitamin D test is a rapid indirect competitive fluorescent immunoassay designed for the quantitative measurement of 25-Hydroxy vitamin D and related hydroxylated metabolites in human serum and plasma (K3EDTA, lithium-heparin and citrate) specimens using the FREND™ System, and the FREND™ AP. Measurements of total 25-hydroxy vitamin D and related hydroxylated metabolites are used to aid in the assessment of vitamin D sufficiency.

The FREND Vitamin D microfluidic flow cartridge is designed for use in the FREND™ System fluorescent immunoassay reader, and the FREND™ AP. The FREND Vitamin D test system is intended for use in clinical laboratories. For in vitro diagnostic use only. The test is not intended for use in point-of-care settings. For professional use only.

2. Summary and explanation of test

Vitamin D is a fat-soluble prohormone known for its role in regulating calcium and phosphorus levels in bone mineralization.^[1,2] Sunlight exposure produces vitamin D via photochemical conversion of 7-dehydrocholesterol in the epidermis and is the primary source of vitamin D.^[3,4] Seasonal changes, amount of exposure, sunscreen use, and skin pigmentation can cause variation in the amount of vitamin D produced in the body. A minor source of

vitamin D can be absorbed from food and vitamin supplements, with an estimated 10-20% absorbed by the body in this manner.^[5,6]

In circulation, 25-OH vitamin D is bound to vitamin D binding protein (VDBP) or albumin at 1000 times higher concentrations than the active form 1,25-(OH)₂-vitamin D.^[7] Additionally, the 25-OH form has a half-life of 2-3 weeks, as compared to the less stable 1,25-(OH)₂ form, which has a half-life of a few hours.^[8] 25-OH vitamin D exists as D2 (ergocalciferol) and D3 (cholecalciferol) isomers, with supplements available for both. Often, total 25-OH vitamin D is measured to assess the sufficiency in a patient and make appropriate clinical decisions.^[9]

3. Principle of the assay

The FRENDS Vitamin D test cartridge is a single-use rapid “competitive” immunoassay utilizing fluorescent nanoparticles in microfluidic flow to capture and quantify total 25-OH vitamin D in serum and plasma specimens. A 35 µL specimen is placed into a Dilution Tube and 70 µL of the diluted sample is placed into the FRENDS Vitamin D pretreatment tube, which contains proprietary reagents. The pretreatment tube is processed in the FRENDS™ AP (Advanced Preparing) device. The AP device mixes and incubates the sample and loads the correct volume onto the FRENDS Vitamin D cartridge. During the incubation labeled 25-OH vitamin D-antibody immune complexes are formed and the AP device transfers a 35 µL sample to the FRENDS Vitamin D cartridge. Vitamin D conjugated fluorescent nanoparticles in the cartridge compete with the 25-OH vitamin D-antibody immune complexes from the pretreatment tube for binding to antibodies immobilized in the test zone. The intensity of fluorescence is inversely proportional to the amount of 25-OH vitamin D in the sample. The FRENDS™ System measures the fluorescence and calculates, displays, stores and optionally prints the 25-OH Vitamin D concentration.

4. Material provided

Qty	Contents	Catalogue number
20	FRENDS Vitamin D cartridge(s)	FRVDAP 020
20	FRENDS Vitamin D pretreatment tube(s)	
20	FRENDS Vitamin D dilution tube(s)	
40	Disposable pipette tip(s)	
01	FRENDS Vitamin D Code chip	
01	FRENDS Vitamin D Package Insert	

One cartridge contains:

Monoclonal mouse anti-vitamin D
25-hydroxyvitamin D
Fluorescent particles

One pretreatment tube contains:

Perfluorohexanoic Acid

One dilution tube contains:

Gold nano-particle conjugation antibody

5. Materials required but not provided

- The FREND™ System
- The FREND™ AP
- Micro-pipette, or any pipette, capable of delivering 35 and 70 µL
- Personal protective equipment and biohazard waste disposal containers

6. Warning and Precautions

Caution: *In the United States, federal law restricts this device to sale by or on the order of a physician.*

- The FREND Vitamin D cartridges are intended for in vitro diagnostic use only.
- Vitamin D cartridges are only to be used on the FREND™ System.
- Vitamin D cartridges are disposable, single use devices. Do not reuse them under any circumstances.
- Allow sealed cartridges to come to room temperature for approximately 15–30 minutes prior to use.
- Cartridges and pretreatment tubes should not be frozen.
- Assure the humidity in the laboratory is in the 10–80% range when tests are run.
- Assure the room temperature remains in the range of 64 °F–77 °F (18–25 °C) when tests are run.
- Avoid cross-contamination between samples by using a new pipette tip for each new specimen.
- Avoid high humidity, direct sunlight or heat in the area used for cartridge storage.
- Inaccurate results are possible if the sample used is contaminated in any way.
- Using specimens containing clotted fibrin could result in erroneous results.

- Over or under loading the cartridge with sample may cause inaccurate results.
- Human specimens are not used in the preparation of this product. However, since human specimens will be used for samples, and since other quality control products in the lab may be derived from human materials please practice Universal Precautions when handling all specimens and controls.
- Do not use the cartridges, pretreatment tubes or dilution tubes beyond the expiration date on the pouch.
- Do not use the cartridge and pretreatment tubes if the pouch is damaged or the seal is broken.
- Perform testing as specified in the Package insert and User manual.
- Keep the cartridge and pretreatment tube sealed in the pouch until ready for use.
- Use the cartridge and pretreatment tube immediately after opening the pouch.
- Wear disposable gloves when handling the cartridges, pretreatment tubes and the samples.
- Wash hands thoroughly and often after handling reagent cartridges or samples.
- Do not ingest the silica gel package found in the cartridge pouch.
- The FREN D Vitamin D assay has been evaluated and no high dose hook effect was observed for Vitamin D concentrations up to 1600 ng/mL.
- Handle specimens in accordance with the OSHA Standard on Bloodborne Pathogens.^[10]

7. Storage and Stability

All unopened materials are stable until the expiration date on the label when stored at the specified temperature. Reagent stability has been demonstrated for twelve months from the date of manufacture.

The expiration date is clearly indicated on the product box and the cartridges

Materials	Catalogue number
Refrigerator temperature storage (2~8 °C)	
FREND Vitamin D cartridges, pouched	FRVDAP 020
FREND Vitamin D pretreatment tubes, pouched	None
FREND Vitamin D dilution tubes	None
Room temperature storage (18~25 °C)	
Pipette tips	None

8. Specimen collection and handling

Human serum and plasma (lithium-heparin, K3EDTA and citrate) samples are suitable for use with FREND Vitamin D cartridges.

Follow instructions detailed in this package insert as well as the specimen collection tube manufacturer's instructions for specimen collection and preparation (including manufacturer's instructions for centrifugation time and speed.)

For serum, a blood sample is collected aseptically without additives by venous puncture. After allowing the sample to clot for 30 minutes at room temperature, the collection tube should be centrifuged for 10 minutes at 3,000 rpm.

For plasma (lithium-heparin, K3EDTA and citrate), a venous blood sample is collected aseptically with the designated additive. After allowing the specimen to sufficiently mix with anticoagulant at room temperature, centrifuge the tube for 10 minutes at 3,000 rpm.

Samples may be stored at 36~46 °F (2~8 °C) for up to 6 hours prior to analysis. If the analysis is scheduled to be done at some later time, the serum or plasma sample should be separated from the red cells and stored frozen at -4 °F (-20°C) or below for future use.

Repeated freeze-thaw cycles should be avoided. Prior to assay, slowly bring frozen samples to room temperature (64~77 °F or 18~25 °C) and mix gently but thoroughly before testing.

For optimal results, avoid grossly hemolytic, lipemic, or turbid specimens. Specimens should be free of aggregated fibrin, red blood cells, or other particulate matter.

9. Procedure

Calibration

There is no need for calibration to be performed by the end user as is generally required on other automated laboratory equipment. All calibration statistics and information have been electronically stored on the FREN D Vitamin D Code chip included in each box of FREN D Vitamin D cartridges. The FREN D Vitamin D Code chip is specific for each lot of FREN D Vitamin D cartridges.

Calibration information should always be checked by running external quality control samples to verify that the results obtained for Vitamin D on the FREN D™ System using the FREN D Vitamin D cartridges meet the laboratory criteria for acceptability.

Code chip installation

Please refer to the FREN D™ System User manual for more detailed instructions relative to the Code chip installation. Abbreviated instructions follow here:

- 1) Insert the FREN D™ System electrical cord into an appropriate outlet.
- 2) Insert the Code chip into the Code chip slot at the rear of the FREN D™ System following the arrows.
- 3) Press the **'Setup'** button on the **'Main'** screen.
- 4) Press the **'Code chip'** button on the **'Setup'** screen.
- 5) The information embedded on the FREN D Vitamin D Code chip is automatically saved on the FREN D™ System.
- 6) When the Code chip installation is completed, press the **'OK'** button to go to the **'Setup'** screen.
- 7) Press the **'Item'** button on the **'Setup'** screen.
- 8) Check the FREN D Vitamin D cartridge lot number and the installation date of the Code chip.
- 9) Press the **'Home'** button to go to the **'Main'** screen to begin running external quality control and patient samples.

Quality control

• FREND™ System QC cartridges

FREND QC Cartridge contains multiple controls to check optic part of the system. By testing QC Cartridge, part of analytical components of the system of (1) laser power, (2) alignment, and (3) mechanical integrity are confirmed. For each day of patient testing, perform QC Cartridge testing. Refer to the quality control procedures section in the User Manual of FREND™ System. In brief, perform QC cartridge testing for the following conditions

- 1) Upon initial setup of the system
- 2) Each day of patient testing,
- 3) When the system has been transported or moved,
- 4) Whenever there is uncertainty about the performance of the system,
- 5) Whenever required by your laboratory's quality control requirements.

• Internal procedural controls

The FREND Vitamin D test cartridge contains a built-in control feature. Fluorescence signal in the reference zone of each cartridge shows: (1) that enough volume is added, (2) that proper flow is obtained, and (3) that the antibody is reactive. If this reference zone signal is missing or lower than threshold, the FREND™ System considers it an incorrect or failed test and produces an error message instead of a test result. In addition, with each cartridge run, the system monitors for (1) flow of sample, (2) speed of sample flow, (3) shelf-life of cartridge components, (4) function of internal barcode scanner, and (5) function of scanner's mechanical components.

• External quality control testing

Commercially available controls from a variety of manufacturers are available that contain 25-OH Vitamin D as a measured analyte. It is recommended that a minimum of two (2) levels of controls be run once per day on days when assaying patient samples on the FREND Vitamin D Test. Or another QC option may be adopted – Individualized Quality Control Plan (IQCP)^[1] is an all-inclusive approach to assuring the quality of the entire testing process. An IQCP includes practices, data and information that each laboratory already uses to ensure quality testing and meet CLIA, beyond testing a certain number of QC materials at a designated frequency.

To ensure that these control procedures are equivalent to CLIA QC regulations and suitable for each laboratory, it is required to establish and document the QC that is appropriate for the test system, testing environment and testing personnel.

Controls should be run according to local, state, federal regulations or accreditation requirements for quality control frequency. If any external quality control sample values are out of the acceptable range, it will be necessary to investigate the problem before reporting patient results to assure there is not an instrument or software malfunction. Do not assay patient samples on the FREND™ System using FREND Vitamin D if quality control results do not fall within the acceptable ranges.

Specimen processing

• Preparation

Remove sufficient FREND Vitamin D cartridges, Sample dilution tubes, and Gold Antibody pretreatment tubes from the refrigerator to test the number of patient samples and required external quality control materials. Allow the tubes and the sealed pouches containing the cartridges to come to room temperature for approximately 15–30 minutes prior to the start of the testing sequence.

When using refrigerated patient samples, remove those from the refrigerator and allow them to come to room temperature prior to testing. If frozen samples will be utilized, be sure these are removed from the freezer, thawed at room temperature and then mixed gently but thoroughly prior to testing. Testing should not begin on previously frozen samples until they have reached room temperature (64–77 °F or 18–25 °C).

There are no other reagents or sample preparations necessary.

• Assay procedure

- 1) Ensure that the FREND Vitamin D cartridge, pretreatment tube, Sample Dilution Tube and specimen are at room temperature (64–77 °F or 18–25°C).
- 2) Open the pouch and place the FREND Vitamin D cartridge into the cartridge tray of the AP device.
- 3) Press **“NEXT”** to close the cartridge tray and open the pretreatment tube tray.
- 4) Transfer 35µL of specimen into the Sample Dilution tube and mix well by vortexing or inverting.
- 5) Transfer a 70 µL aliquot of the diluted specimen to the FREND Vitamin D pretreatment tube.

Caution: *Once the sample is added to the pretreatment tube. Do not invert the tube.*

- 6) Insert the pretreatment tube into the tube hole in the FREND AP pretreatment tube tray. Refer to the FREND AP User manual for

- complete operating instructions.
- 7) Press the **"NEXT"** button. The pretreatment tray will close and the first incubation step (15 minutes) will begin.
 - 8) After the first incubation is complete, 35 µL of mixed sample will be loaded onto the cartridge and the second incubation step (2 minutes) will begin.
 - 9) When both incubation steps are completed, the cartridge tray will open and the cartridge will be ready to be inserted into the FREND™ System.
 - 10) Press the **'Test'** button on the **'Main'** screen of the FREND™ System.
 - 11) The system moves to the Patient ID screen automatically.
 - 12) Type the Patient ID and press the **'Enter'** button to begin the test.
 - 13) Insert the cartridge into the cartridge slot using the cartridge arrow as a guide.

Caution: Check the direction of the cartridge before insertion and assure the insertion is complete.

- 14) When the reaction in the cartridge is completed, the FREND™ System will automatically begin the reading process.
- 15) When the measurements are completed, the cartridge will automatically be expelled and the results displayed.

Caution: Do not disconnect power cord or shut off power on the FREND™ System while a cartridge is in the reading chamber. This may cause a system error.

- 16) If the FREND™ System is connected to the optional printer, press the **'Print'** button and the results will be output on the printer paper.
- 17) For more detailed instructions, please refer to the FREND™ System User manual.

10. Procedural notes

If the Vitamin D concentration of a specimen is found to be greater than the linearity limit of the assay of 110.00 ng/mL, report any value above 110.00 ng/mL as "greater than 110.00 ng/mL". There is no specific diluent provided in the FREND Vitamin D cartridge.

11. Calculation of results

The FREND™ System performs all sample and reagent handling operations automatically within the cartridge once the sample has been loaded into the sample inlet in the cartridge and the cartridge placed into the FREND™ System. The rate of fluorescence produced by the reaction is read at various intervals during the analysis process, blank readings are subtracted after

which the net rate is automatically converted to Vitamin D concentration in ng/mL based upon information stored on the FREND Vitamin D Code chip. This result is then output on the screen and to the optional printer. It is also stored in memory on the FREND™ System.

Screen displayed for various concentration scenarios

Displayed result	Description
	Vitamin D concentration Less than 10.00 ng/mL
	Vitamin D concentration Not less than 10.00 ng/mL and not higher than 110.00 ng/mL
	Vitamin D concentration Higher than 110.00 ng/mL

12. Limitations of the procedure

- 1) When used for diagnostic purposes, the results obtained from this assay should be used in conjunction with other data (e.g. symptoms, results of other tests, clinical impressions, medical history, therapy, etc.)
- 2) The FREND™ System paired with a FREND Vitamin D cartridge, is programmed to report 110.00 ng/mL as the highest concentration of Vitamin D measurable. The lowest measurable concentration is 13.0 ng/mL, the assay limit of detection.
- 3) Specimens from patients with heterophilic antibodies, such as anti-mouse (HAMA), anti-goat (HAGA), or anti-rabbit (HARA) antibodies, may show falsely elevated or depressed values or may result in an incomplete test.^[12, 13] Patients routinely exposed to animals or animal serum products can be prone to these types of heterophilic interferences.
- 4) Patients taking the drugs containing Paricalcitol (e.g., Zemplar) should not be tested by this assay.
- 5) Certain medications may interfere with assay performance. All results should be interpreted with respect to the clinical picture of the patient.
- 6) Although hemolysis has an insignificant effect on the assay, hemolyzed samples may indicate mistreatment of a specimen prior to assay and results should be interpreted with caution.
- 7) Lipemia has an insignificant effect on the assay except in the case of gross lipemia where interference with the lateral flow of the sample in the cartridge may occur.
- 8) The concentration of Vitamin D in a given sample determined with assays from different manufacturers can vary due to differences in assay methods, calibration, and antibody specificity.
- 9) Please refer to the Specimen Collection and Handling, Warnings and Precautions, Storage and Stability, and Procedural Notes sections in this Package insert.
- 10) FREND Vitamin D has not been evaluated in point-of-care settings.
- 11) FREND Vitamin D is to be used in licensed clinical laboratories with trained technologists.

13. Performance characteristics

Precision

Clinical and Laboratory Standards Institute (CLSI) document EP5-A2 was utilized as a guideline for precision studies on FRENDD Vitamin D. Standard material with three different concentration was used for precision test. One operator tested with one lot for 20 days. One operator tested with three different lots for 4 times. Three different operators tested with same lot for 4 times.

1) One operator with one lot for 20 days

Standard material	Conc. (ng/mL)	Within-run		Between-run		Between-day		Within-laboratory	
		SD	CV(%)	SD	CV(%)	SD	CV(%)	SD	CV(%)
Low	30	2.458	8.1	0.895	2.9	0.689	2.3	2.705	8.9
Medium	60	0.146	4.9	0.040	1.3	0.021	0.7	0.153	5.1
High	100	4.865	4.9	0.801	0.8	1.336	1.3	2.108	5.1

2) One operator with three different lots for 4 times


Standard material	Conc. (ng/mL)	Lot. 503904	Lot. 503905	Lot. 503906	Mean	SD	CV(%)
		Mean	Mean	Mean			
Low	30	31.14	31.913	30.553	31.292	0.682	2.2
Medium	60	58.98	59.473	60.68	59.711	0.875	1.5
High	100	95.528	100.645	98.085	98.086	2.559	2.6

3) Three different operators with one lot for 4 times

Standard material	Conc. (ng/mL)	Operator1	Operator2	Operator3	Mean	SD	CV(%)
		Mean	Mean	Mean			
Low	30	29.430	28.890	30.755	29.692	0.960	3.2
Medium	60	55.775	56.025	60.923	57.574	2.902	5.0
High	100	98.933	101.598	99.320	99.950	1.440	1.4

Dilution Linearity

The range of linearity in serum was established using CLSI document EP6-A. FREN D Vitamin D was found to have linearity within the reportable range of 8.00 ng/mL ~ 130.00 ng/mL with a mean recovery of 100±10%. The dilution linearity study was performed by diluting a high concentration Vitamin D specimen with a low concentration Vitamin D specimen.

Analytical sensitivity


The limit of blank (LoB) and limit of detection (LoD) was determined using guidelines found in CLSI document EP17-A. LoB was determined from 60 replicate measurements using a calibrator A (vitamin D depleted serum). LoD was determined using 12 replicates measurements of five low patient samples.

LoB	LoD	LoQ
3.85 ng/mL	6.53 ng/mL	10.00 ng/mL


Comparative analysis

FREND Vitamin D was compared to the predicate device using guidelines outlined in CLSI document EP09-A2-IR. Samples (n=120) were measured in duplicate on both systems. Linear regression analysis demonstrated a correlation (r) of 0.9811 and 0.9774.

1) FREND Vitamin D VS Predicate device A

2) FREND Vitamin D VS Predicate device B

Cross-reactivity

The following substances were evaluated for potential cross-reactivity with the FRENED Vitamin D at three concentrations. Testing was done according to the CLSI protocol EP7-A2. No significant cross-reactivity was found.

No.	Potential cross reactants	Concentration
1	Vitamin D2(Ergocalciferol)	1,000 ng/mL
2	Vitamin D3(Cholecalciferol)	1,000 ng/mL
3	1,25-(OH) ₂ -Vitamin D2	100 ng/mL
4	1,25-(OH) ₂ -Vitamin D3	100 ng/mL
5	3-epi-25-(OH) Vitamin D3	100 ng/mL
6	25-(OH) Vitamin D2	25 ng/mL
7	25-(OH) Vitamin D3	25 ng/mL

Interference

The interference study was performed as recommended in the CLSI EP7-A2 protocol using three concentrations of Vit D. Recovery within 90% to 110% of the expected Vit D was considered as lack of interference. No interference by the substances below was found.

No.	Potential interfering	Concentration
1	Bilirubin	20.0 mg/dL
2	Biotin	30.0 ng/dL
3	Triglycerides	370.0 mg/dL
4	Human albumin	12.0 g/dL
5	Hemoglobin	200.0 mg/dL

Matrix effect

The heparin-plasma, EDTA-plasma, and Citrate-plasma can be used for the FRENED Vitamin D assay. The correlation values are $R^2=0.9816$ in serum with heparin-plasma, $R^2=0.9739$ in serum with EDTA-plasma, and $R^2=0.9797$ in serum with Citrate-plasma in each test and the values shows the high correlation.

14. References


- 1) Holick, MF. High prevalence of vitamin D inadequacy and implications for health. *Mayo Clinic Proceedings*. 2006, 81 (3): 353-373.
- 2) Bolland, MJ et al. The effect of vitamin D supplementation on skeletal, vascular, and cancer outcomes: a trial sequential meta-analysis. *The Lancet Diabetes & Endocrinology*. 2014, 2 (4): 307-320.
- 3) Wolf, G. The discovery of vitamin D: the contribution of Adolf Windaus. *Journal of Nutrition*. 2004, 134 (6): 1299-1302.
- 4) Holick, MF. Vitamin D: importance in the prevention of cancers, type 1 diabetes, heart disease, and osteoporosis. *The American Journal of Clinical Nutrition*, 2004, 79 (3): 362-371.
- 5) Calvo, MS et al. Vitamin D intake: a global perspective of current status. *Journal of Nutrition*, 2005, 135 (2): 310-316.
- 6) Ross, AC et al. Dietary Reference Intakes for Calcium and Vitamin D. Washington DC: *National Academies Press* 2011, p. 435.
- 7) Svasti, J et al. Molecular basis for the three major forms of human serum vitamin D binding protein. *Biochemistry* 1979, 18 (8): 1611-1617.
- 8) Holick, MF et al. Vitamin D2 is as effective as vitamin D3 in maintaining circulating concentrations of 25- hydroxyvitamin D. *The Journal of Clinical Endocrinology and Metabolism*, 2008, 93 (3): 677-681.
- 9) Maxmen, A. Nutrition advice: the vitamin D-lemma. *Nature* 2011, 475 (7354): 23-25.
- 10) Chao, E.L.; Henshaw, J.L. Occupational Safety and Health Administration: Model Plans and Programs for the OSHA Bloodborne Pathogens and Hazard Communications Standards. *OSHA 3186-06R*, 2003.
- 11) Centers for Medicare & Medicaid Services. CLIA Individualized Quality Control Plan – Considerations when deciding to develop an IQCP. *CLIA Brochure #12*. November 2014.
- 12) Schroff, R.W.; Foon, K.A.; Beatty, S.M.; Oldham, R.K.; Morgan, A.C. Human Anti-Murine Immunoglobulin Responses in Patients Receiving Monoclonal Antibody Therapy. *Cancer Research*. 1985, 45: 879-885.
- 13) Boscatto, L.M.; Stuart, M.C. Heterophilic Antibodies: A Problem for All Immunoassays. *Clinical Chemistry*, 1988, 34 (1): 27-33.

Glossary of Symbols

	Caution, warning, Consult accompanying documents
	Catalogue number/Reference number
 www.nanoentek.com/eifu.php	Consult Instructions for Use An electronic instructions for use (eIFU) indicator (website address) may accompany the symbol when used to indicate an instruction to consult an eIFU.
	Lot number/Batch number
	Use by YYYY-MM-DD or YYYY-MM
	Manufacturer
	CE marking
	<i>In vitro</i> diagnostic medical device
	Temperature limitation
	Contains sufficient for <n> tests
	Do not reuse
	Do not use if package is damaged
	For prescription use only CAUTION: Federal (U.S.) law restricts this device to sale by or on order of a physician.
	US Corporation
	Patient ID
	Result
	Sample Drop
	Authorized representative in the European Community

UK Representative	Authorized representative in United Kingdom
CH REP	Authorized representative in Switzerland
BRH	Authorized representative in Brazil

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14. Referenzen

- 1) Holick, M.F. Hohe Prävalenz des Vitamin-D-Mangels und Auswirkungen auf die Gesundheit. *Mayo Clinic Proceedings*. 2006, 81 (3): 353-373.
- 2) Bolland, M.J., et al. The effect of vitamin D supplementation on skeletal, vascular, and cancer outcomes: a trial sequential meta-analysis. *The Lancet Diabetes & Endocrinology*. 2014, 2 (4): 307-320.
- 3) Wolf, G. Die Entdeckung von Vitamin D: der Beitrag von Adolf Windaus. *Journal of Nutrition*. 2004, 134 (6): 1299-1302.
- 4) Holick, M.F. Vitamin D: Die Bedeutung für die Prävention von Krebs, Typ-1-Diabetes, Herzerkrankungen und Osteoporose. *The American Journal of Clinical Nutrition*, 2004, 79 (3): 362-371.
- 5) Calvo, MS et al. Vitamin D intake: a global perspective of current status. *Journal of Nutrition*, 2005, 135 (2): 310-316.
- 6) Ross, AC et al. Dietary Reference Intakes for Calcium and Vitamin D. Washington DC: *National Academies Press* 2011, p. 435.
- 7) Svasti, J et al. Molecular basis for the three major forms of human serum vitamin D binding protein. *Biochemistry* 1979, 18 (8): 1611-1617.
- 8) Holick, MF et al. Vitamin D2 is as effective as vitamin D3 in maintaining circulating concentrations of 25- hydroxyvitamin D. *The Journal of Clinical Endocrinology and Metabolism*, 2008, 93 (3): 677-681.
- 9) Maxmen, A. Nutrition advice: the vitamin D-lemma. *Nature* 2011, 475 (7354): 23-25.
- 10) Chao, E.L.; Henshaw, J.L. Occupational Safety and Health Administration: Model Plans and Programs for the OSHA Bloodborne Pathogens and Hazard Communications Standards. *OSHA 3186-06R*, 2003.
- 11) Centers for Medicare & Medicaid Services. CLIA Individualized Quality Control Plan – Considerations when deciding to develop an IQCP. *CLIA Brochure #12*. November 2014.
- 12) Schroff, R.W.; Foon, K.A.; Beatty, S.M.; Oldham, R.K.; Morgan, A.C. Human Anti-Murine Immunoglobulin Responses in Patients Receiving Monoclonal Antibody Therapy. *Cancer Research*. 1985, 45: 879-885.
- 13) Boscatto, L.M.; Stuart, M.C. Heterophilic Antibodies: A Problem for All Immunoassays. *Clinical Chemistry*, 1988, 34 (1): 27-33.

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14. Riferimenti

- 1) Holick, MF. High prevalence of vitamin D inadequacy and implications for health. *Atti della Mayo Clinic*. 2006, 81 (3): 353-373.
- 2) Bolland, MJ et al. The effect of vitamin D supplementation on skeletal, vascular, and cancer outcomes: a trial sequential meta-analysis. *The Lancet Diabetes & Endocrinology*. 2014, 2 (4): 307-320.
- 3) Wolf, G. The discovery of vitamin D: the contribution of Adolf Windaus. *Journal of Nutrition*. 2004, 134 (6): 1299-1302.
- 4) Holick, MF. Vitamin D: importance in the prevention of cancers, type 1 diabetes, heart disease, and osteoporosis. *American Journal of Clinical Nutrition*, 2004, 79 (3): 362-371.
- 5) Calvo, MS et al. Vitamin D intake: a global perspective of current status. *Journal of Nutrition*, 2005, 135 (2): 310-316.
- 6) Ross, AC et al. Dietary Reference Intakes for Calcium and Vitamin D. Washington DC: *National Academies Press* 2011, p. 435.
- 7) Svasti, J et al. Molecular basis for the three major forms of human serum vitamin D binding protein. *Biochemistry* 1979, 18 (8): 1611-1617.
- 8) Holick, MF et al. Vitamin D2 is as effective as vitamin D3 in maintaining circulating concentrations of 25- hydroxyvitamin D. *The Journal of Clinical Endocrinology and Metabolism*, 2008, 93 (3): 677-681.
- 9) Maxmen, A. Nutrition advice: the vitamin D-lemma. *Nature* 2011, 475 (7354): 23-25.
- 10) Chao, E.L.; Henshaw, J.L. Occupational Safety and Health Administration: Model Plans and Programs for the OSHA Bloodborne Pathogens and Hazard Communications Standards. *OSHA 3186-06R*, 2003.
- 11) Centers for Medicare & Medicaid Services. CLIA Individualized Quality Control Plan – Considerations when deciding to develop an IQCP. *CLIA Brochure #12*. November 2014.
- 12) Schroff, R.W.; Foon, K.A.; Beatty, S.M.; Oldham, R.K.; Morgan, A.C. Human Anti-Murine Immunoglobulin Responses in Patients Receiving Monoclonal Antibody Therapy. *Cancer Research*. 1985, 45: 879-885.
- 13) Boscatto, L.M.; Stuart, M.C. Heterophilic Antibodies: A Problem for All Immunoassays. *Clinical Chemistry*, 1988, 34 (1): 27-33.

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NanoEntek, Inc.
MT Promedt Consulting GmbH